

# Detection of PD-L1 and PD-L2 on Circulating Tumor Cells (CTCs) Using Microfluidic Based Chipcytometry

Jinkai Teo<sup>1</sup>, Meihui Tan<sup>2</sup>, Anja Mirenska<sup>3</sup>, Janice Oh<sup>1</sup>, Lewis Hong<sup>1</sup>, Richard Wnek<sup>4</sup>, Jan Detmers<sup>3</sup>, Ali Asgar S. Bhagat<sup>2</sup>, Chih-Liang Chin<sup>1</sup>, David Skibinski<sup>1</sup>

<sup>1</sup>Translational Biomarkers, Translational Medicine Research Centre, Merck Research Laboratories, MSD, Singapore

<sup>2</sup>Clearbridge BioMedics Pte Ltd, 81 Science Park Drive, The Chadwick, #02-03, Singapore Science Park 1 Singapore 118257

<sup>3</sup>Clinical Biomarkers, Zellkraftwerk, Bosestrasse 4, D-04109 Leipzig, Germany


<sup>4</sup>Translational Molecular Biomarkers, Merck Research Laboratories, Merck & Co. Inc., Rahway NJ 07065, USA



## Abstract

- Recent advances in cancer therapy have demonstrated the potential of the immune system in cancer control and rejection.
- Prominent amongst these approaches has been the success of anti-PD-1 immunotherapy, to break the strong inhibitory signal, transmitted by tumor specific ligands such as PDL-1, to the PD-1 immune modulatory receptor expressed on T-cells.
- PD-L1 expression on the tumor is a clinically validated biomarker of therapeutic response to anti-PD-1 immunotherapy. However, obtaining tumor biopsies for PD-L1 interrogation is an invasive procedure not suited for frequent longitudinal monitoring during cancer therapy. Furthermore, tumor heterogeneity for PD-L1 expression may not accurately capture the PD-L1 status of the whole tumor burden in a single biopsy.
- An alternative, minimally invasive, approach is the analysis of blood samples for circulating tumor cells (CTCs) which have broken away from the tumor and entered the periphery.
- We describe the development of an assay workflow to detect and characterize circulating tumor cells in peripheral blood samples. Our approach uses a sized-based microfluidic enrichment technique, and subsequent characterization with microfluidic based cytometry (Chipcytometry).


## Workflow



## CTC enrichment by ClearCell® FX System



- Microfluidic chip enriches CTCs based on size
- Label free isolation
- Process blood volume of 7.5 mL
- Processing in approx. 1 hour




- Smaller cells (RBCs ~8 μm and leucocytes ~8–15 μm) are affected by the Dean drag and migrate to the outer wall
- Larger CTCs (~15–20 μm) experience strong inertial lift forces and are focused along the inner wall


## Specificity of anti-human PD-L1 and PD-L2 monoclonal antibodies determined by Flow Cytometry

- Specificity of monoclonal antibodies against PD-L1 and PD-L2 was determined on non-small cell lung cancer cell lines, A549 and H1975, using flow cytometry.
- A549 cells were incubated with IFN-γ for 24 hours prior to staining to upregulate PD-L1 and PD-L2. Cells were fixed prior to staining.
- The data demonstrates the presence of PD-L1 on A549 and H1975 cells, and PD-L2 on A549 cells.
- The specificity of the staining was demonstrated by the fact that no signal was detected for the isotype control on both A549 and H1975 cells.


Flow cytometry histograms:



## CTC characterization by Chipcytometry



- Cells immobilized on microfluidic chip
- Fixation of immobilized cells and long term storage at 4°C



- Phenotyping by staining with monoclonal antibodies conjugated to a fluorophore and detection using CYTOBOT system
- Sequential staining of each antibody permits unlimited set of markers



- Individual cells are identified and based on their position
- Multidimensional information is layered to create a profile for each cell


Marker	Relevance
DAPI	Validated marker of CTC identification*
CD45	Validated marker of CTC identification*
Cytokeratin	Validated marker of CTC identification*
EpCAM	Validated marker of CTC identification*
CD3	Negative on CTCs
CD15	Negative on CTCs
PD-L1	Immune checkpoint regulator
PD-L2	Immune checkpoint regulator
Vimentin	Marker of Epithelial – mesenchymal transition

\* FDA validated CellSearch system identify CTCs as DAPI+ / CD45- / Cytokeratin+ / EpCAM+


## Detection of PD-L1 and PD-L2 expression by Chipcytometry

- A549, H1975 and whole blood was analyzed for PD-L1 and PD-L2 expression using the CYTOBOT chipcytometry system. A549 cells were incubated with IFN-γ for 48 hours prior to staining to upregulate PD-L1 and PD-L2. Histograms translated from images.
- Individual A549 and H1975 cells positive for PD-L1, and individual A549 cells positive for PD-L2, could be observed.


### Detection of PD-L1



### Detection of PD-L2



Histograms derived from chipcytometry image fluorescence:



## Conclusions

- Assay workflow to detect and characterise CTCs
- Specificity of PD-L1 and PD-L2 antibodies demonstrated by flow cytometry
- Established feasibility of PD-L1 and PD-L2 detection by chipcytometry
- Evaluation of workflow in currently on-going experiments with spiked samples and samples from breast cancer patients